

Ingredients for Mobility

Below are the parts that are included when a wheelchair is shipped. Using the Assembly Manual and included tools, a chair can be easily assembled anywhere in the world.

1



Plastic Lawn Chair (manufactured specifically to meet our needs)

Shipping box for wheels and other hardware

Footrest for the front of the wheelchair

Nuts and Bolts

Steel Frame and Assembly Tools

6



Completed Wheelchair

2



Castors

Brakes

Bicycle Pump

Patch Kit

Axle for main wheels

Axle Supports

24 inch Mountain Bike Tires

5



Wheelchairs are assembled on-site

3




Footplate

Harness

Pressure Cushion

Assembly Instructions

4



Wheelchairs arrive unassembled