



Ingredients for Mobility

Below are the parts that are included when a wheelchair is shipped. Using the Assembly Manual and included tools, a chair can be easily assembled anywhere in the world.

1



Plastic Lawn Chair
(manufactured specifically
to meet our needs)




Shipping box for wheels
and other hardware



Footrest for the front
of the wheelchair



Nuts and Bolts



Steel Frame and
Assembly Tools

6




Completed Wheelchair

2



Castors



Brakes




Bicycle Pump




Patch Kit



Axle for
main wheels



Axle
Supports



24 inch Mountain Bike Tires

5




Wheelchairs are assembled on-site


3



Footplate



Harness



Pressure Cushion



Assembly
Instructions

4



Wheelchairs arrive unassembled