

# OUR WHEELCHAIRS

Learn  
More



For millions of people living with disabilities in the developing world, a wheelchair can open the door to independence, opportunity, and full participation in community life. Cost-effective, durable models by Free Wheelchair Mission help make that possible.

Since 2001, Free Wheelchair Mission has provided 1.5 million wheelchairs to people living with disabilities across 95 countries. Provided at no cost to recipients through the gifts of generous supporters, our wheelchairs provide renewed dignity, independence, freedom, and joy.

## GEN 4 JOY™

Introduced 2025



## GEN 3

Introduced 2013



## GEN 2

Introduced 2009



## GEN 1

Introduced 2001  
Phased Out 2020



### **GEN 4 JOY: More Sizes. Expanded Features. BRINGS JOY.**

Our newest wheelchair brings more sizes, from XS to XXL. Packed with new features for greater comfort, improved posture, smoother maneuverability, and long-lasting durability. A lower seat with adjustable depth and customizable footrest angles enhance support, while its foldable, compact frame improves transport. It introduces fast, easy assembly and adjustments.

### **GEN 3: Foldable. Transport-Friendly. EXPANDS FREEDOM.**

Introduced foldability for easier storage, transport, and use in tight spaces.

### **GEN 2: Multiple Sizes. Adjustable Fit. INCREASES INDEPENDENCE.**

Brought greater adjustability for a personalized fit for smaller recipients or those with postural support needs.

### **GEN 1: Durable. Reliable. Rugged. RENEWED DIGNITY.**

Our original model, built for durability, provided reliable mobility in rugged terrain. GEN 1 came in a single size and was phased out in 2020.

[FreeWheelchairMission.org](https://FreeWheelchairMission.org)

15279 Alton Parkway, Suite 300, Irvine, CA 92618 • 949.273.8470 • [info@freewheelchairmission.org](mailto:info@freewheelchairmission.org)



# OUR WHEELCHAIRS

Learn  
More



All four of our wheelchair models are designed to be durable, withstanding the daily challenges faced in parts of the developing world. Our thoughtful approach to engineering and distribution allows us to maximize impact worldwide—providing mobility along with **renewed dignity, independence, freedom, and joy.**

	GEN 1 Phased Out 2020	GEN 2	GEN 3	GEN 4 JOY
Introduced	2001	2009	2013	2025
Sizes	one size	S, M, L, XL	S, M, L, XL	XS, S, M, L, XL, XXL
Adjustable & Removable Footrests		✓	✓	✓
Foldable			✓	✓
Quick Release Wheels				✓
Upright Seat				✓
Wide Armrests				✓
Overall Height	25.1"	31.4"-37.3"	29.8"-37.4"	33.25"-38.5"
Overall Depth	40.5"	43.5"	43.6"	37.4"-45.7"
Seat Width	17.3"	13.5"-20"	13.5"-19"	11"-21.5"
Folded Width	N/A	N/A	15"	13"
Caster Size	8"	8"	8"	6"-8"
Rear Wheel Size	24"	26"	26"	26"
User Weight Capacity	220 lbs	220 lbs	220 lbs	165-275 lbs
Wheelchair Weight	38.4 lbs	36.4 lbs	45.6 lbs	42-48.5 lbs
Wheelchairs Per Container	520	510 Average		up to 550
Estimated Useful Life		Baseline		up to 2x longer!

**To learn more or inquire about supporting our mission, scheduling media interviews, or becoming a distribution partner, call or email us today!**

[FreeWheelchairMission.org](http://FreeWheelchairMission.org)

15279 Alton Parkway, Suite 300, Irvine, CA 92618 • 949.273.8470 • [info@freewheelchairmission.org](mailto:info@freewheelchairmission.org)

